



YP·DENTAL

# Dental Glass Ceramic Blocks



**Anyang Yingpai Dental Material Co., LTD**

Web: <https://ypdental.com>

E-mail: [business@ypdental.com](mailto:business@ypdental.com)

Add: Shangguan Town, Hua County, Henan province,  
China (Mainland)

*Anyang Yingpai Dental Material Co., LTD*

## Enables efficient CAD/CAM processes and precise restoration


Dental glass ceramic blocks select high-quality raw material, with fine craftsmanship, fast and simple crystallization process, excellent bending strength and processability as well as excellent aesthetic restoration effect, becoming the first choice for CAD/CAM chairside material.



- Simple and fast process, convenient treatment for patients
- Rapid crystallization, the best choice for immediate restoration.
- High strength, suitable for various cases.
- High transparency, matching the natural color of the teeth.



## Transparency



## Specification Classification

	<b>C14</b>	Size(mm)	pcs/box		<b>B40</b>	Size(mm)	pcs/box
		18*15*13	5pcs			40*15*15	4pcs

## Types

Colors	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	BL1	BL2	BL3	BL4	OM1	OM2	OM3
HT (C14)	•	•	•	•	•	•	•			•	•						•	•	•	•	•	•	•
LT (C14)	•	•	•	•	•	•	•			•	•			•	•		•	•	•	•	•	•	•
HT (B40)	•	•	•	•	•	•	•			•	•						•	•	•	•	•	•	•
LT (B40)	•	•	•	•	•	•	•			•	•			•	•		•	•	•	•	•	•	•


### Homogeneity

- Smooth surface
- Uniform color
- No impurities and foreign matter




### Bending strength

The average biaxial bending strength of crystallized YPDENTAL glass ceramics is 637Mpa, which is higher than the standard strength value set by ISO 6872:2008 dental ceramic material. Its superior strength can be used to make three-unit bridge for the molar teeth.



### Recommended Indication



Veneer   Inlay   Onlay   Anterior crown   Posterior crown   Monolithic Bridge (≤3 Units)

### Physical properties

Density(g/cm <sup>3</sup> )	2.4-2.7
Vickers hardness(MPa)	5400±400
Fracture toughness(MPa·m <sup>1/2</sup> )	2.55
Coefficient of thermal expansion(10 <sup>-6</sup> K <sup>-1</sup> )	9.8±0.5
Crystallization temperature(°C)	840-850
Chemical solubility(μg·cm <sup>-2</sup> )	30.3
Radioactivity(Bq/g)	< 0.019

### Chemical composition

SiO <sub>2</sub>	61%-71%
Li <sub>2</sub> O	11%-17%
ZrO <sub>2</sub>	1%-5%
Other oxides	6%-30%

### Biological properties

Cytotoxicity test	Cytotoxicity level 0
Short-term systemic toxicity test (oral route)	no systemic toxicity
Sensitization test	no sensitization
Hemolysis test	hemolysis rate <5%
Ames test	negative for mutagenesis
Oral mucosa irritation test	no oral mucosa irritation
Subchronic Systemic Toxicity Test	not cause subchronic systemic toxicity

### Crystallization curve

One-stage crystallization procedure	Starting temperature	Drying time	Heating rate	Holding temperature	Residence time	Furnace temperature	Vacuum start	Vacuum end
	400°C	30s	40°C/min	840°C	7min	700°C	550°C	840°C
Two-stage crystallization process (Recommended for use)	Starting temperature	Drying time	Heating rate	Holding temperature	Residence time	Heating rate	Holding temperature	Residence time
	400°C	6min	60°C/min	820°C	10s	30°C/min	840°C	7min
	Furnace temperature	Vacuum start	Vacuum end					
	700°C	550/820°C	820/840°C					