

YP Dental W500 Milling Machine

Shipping List

No.	Product Name	Specification	Qty	Picture
1	W500 milling machine	W500	1	
2	Power Cord	European Standard_220V	1	
3	Spindle Water Cooler		1	
4	Water Pump		1	
5	Power Box Connector	(Spare)	1	
6	Cutting Fluid Tank		1	

7	Glass Ceramic Fixture		1	
8	Pre-mill Abutment Fixture		1	
9	Tool Wrench		1	
10	Spindle Maintenance Tool		1	
11	Tip Set Screw for Glass Ceramic	M4*8	20	
12	Tip Set Screw for Glass Ceramic	M4*12 (Spare)	20	
13	Screw for Fixture	M5*16	40	
14	Screw for Pre-mill	M5*12	40	

15	Screw Washers for Pre-mill	M5	20	
16	Screw for Disk Fixture	M4*14	20	
17	Screw for Cutting Fluid Filter	M8*10	4	
18	Cutting Fluid Valve	$\Phi 10$	1	
19	Air Shut-off Valve	$\Phi 8$	1	
20	Tee	$\Phi 8$ to $\Phi 8$	1	
21	Tee	$\Phi 8$ to $\Phi 6$	1	
22	Tee	$\Phi 10$ to $\Phi 8$	1	

23	Two-way Valve	$\Phi 10$ to $\Phi 8$	1	
24	Air Filter		1	
25	Cutting Fluid Filter		1	
26	Vacuum Hose		1	
27	Air Hose	$\Phi 8$ (Transparent)	2m	
28	Coolant Hose	$\Phi 8$ (Red)	2m	
29	Cutting Fluid Inlet Hose	$\Phi 10$ (Red)	4m	
30	Cutting Fluid Drain Hose		1	

31	Brush		1	
32	Test Wax Disk	Φ98x10mm	2	
33	Fuse		5	
34	CAM Software USB Flash Drive		1	
35	Screwdriver		1	
36	Pliers		1	
37	Milling Bur for Metal_T1	R1.5*6D*50L	1	
38	Milling Bur for Metal_T2	R1.0*6D*50L	1	

39	Milling Bur for Metal_T3	R0.5*6D*50L	1	
40	Milling Bur for Metal_T4	R1.0*6D*50L	1	
41	Milling Bur for Metal_T5	R0.75*6D*50L	1	
42	Milling Bur for Metal_T6	R1.0*6D*50L	1	
43	Milling Bur for Metal_T7	R0.5*6D*50L	1	
44	Milling Bur for Metal_T8	R0.75*6D*50L	1	
45	Milling Bur for Metal_T9	R0.75*6D*50L	1	
46	Milling Bur for Metal_T10	R1.1*6D*50L	1	

47	Milling Burs for Glass Ceramic_T12	R1.25*6D*50L	1	
48	Milling Burs for Glass Ceramic_T13	R0.5*6D*50L	1	
49	Milling Burs for Glass Ceramic_T14	R0.3*6D*50L	1	
50	Milling Burs for PMMA_T15	R0.3*6D*50L	1	
51	Milling Burs for Zirconia_T1	R1.0*6D*50L	1	
52	Milling Burs for Zirconia_T2	R0.5*6D*50L	1	
53	Milling Burs for Zirconia_T3	R0.3*6D*50L	1	
54	Milling Burs for Zirconia_T4	R0.15*6D*50L	1	

55	Milling Burs for Zirconia_T5	R0.75*6D*50L	1	
----	---------------------------------	--------------	---	---